SYSTEM PACKAGES

Everything you need to get started!



Inlight Therapy polychromatic light therapy devices are FDA CLEARED for Increased Circulation, and Pain Relief.



This package includes three (3) polychromatic therapy pads providing the power of 500 polychromatic LEDs.

- 6 PORT PRO Controller
- FACEMASK 104 Blue & Red PLT
- LOCAL 132 Red & NIR PLT
- BODY 264 Red & NIR PLT



Inlight Therapy polychromatic light therapy devices are FDA CLEARED for Increased Circulation, and Pain Relief.





Inlight Therapy polychromatic light therapy devices are FDA CLEARED for Increased Circulation, and Pain Relief.

2 PORT | 2 PAD PAIN

Experience the usability of the 2/PORT 2 Pad Pain portable system. Two popular pads are easy to use and allow multiple pad placement. Get Polychromatic Light Therapy (PLT) on the go with Inlight Progressive Multi-Pulse[™] technology.



Inlight Therapy polychromatic light therapy devices are FDA CLEARED for Increased Circulation, and Pain Relief.

LED PADS nlight Interchangeable between all Inlight controllers. Controller required and sold separately. 🔵 BLUE | 470 nm 🛛 😑 RED | 630 nm 🛛 🔿 NIR | 850 nm AVAILABLE ONLY: RED/NIR AVAILABLE ONLY: RED/NIR PODIATRIC BOOT PODIATRIC BOOT SINGLE 122 PAIR 122 (x2) Neoprene flexibility for ease of application Neoprene flexibility for ease of application anywhere on the body and on the foot. anywhere on the body and on the foot. · Comfortably secured by a Velcro strap. • Comfortably secured by a Velcro strap. • 122 diodes of the polychromatic blend of red • 122 diodes of the polychromatic blend and infrared photon energy in each pad. of red and infrared photon energy. • 52 Red diodes (630 nm) in each pad. • 52 Red diodes (630 nm). • 70 NIR diodes (850 nm) in each pad. • 70 NIR diodes (850 nm). 🛑 RED | 🔵 NIR 🛑 RED | 🔵 NIR OPTIONS: BLUE/NIR, RED/NIR, or BLUE/RED/NIR AVAILABLE ONLY: RED/NIR **F-PAD 263 BODY 264** • 8"x15" neoprene flexibility for ease of • Flexible "T" shape neoprene therapy pad for ease of application anywhere on the body. application anywhere on the body. • Comfortably secured by a Velcro strap. Comfortably secured by a Velcro strap. • 264 diodes of the polychromatic blend • 263 diodes of the polychromatic blend of red and infrared photon energy. of red and infrared photon energy. • 120 Red diodes (630 nm). • 120 Red diodes (630 nm). • 144 NIR diodes (850 nm). • 143 NIR diodes (850 nm). 🛑 RED | 🔵 NIR 🔵 BLUE | 🛑 RED | 🔵 NIR OPTIONS: BLUE/NIR, RED/NIR, or BLUE/RED/NIR OPTIONS: BLUE/NIR, RED/NIR, or BLUE/RED/NIR PAINBUSTER II 180 **LOCAL 132** 5"x16" neoprene flexibility for ease of application 8"x9" neoprene flexibility for ease of anywhere on the body. application anywhere on the body. Comfortably secured by a Velcro strap. Comfortably secured by a Velcro strap. • 180 diodes of the polychromatic blend of red • 132 diodes of the polychromatic blend and infrared photon energy. of red and infrared photon energy. 80 Red diodes (630 nm) in each pad. • 60 Red diodes (630 nm). • 72 NIR diodes (850 nm). • 100 NIR diodes (850 nm) in each pad. 🔵 BLUE | 🛑 RED | 🔾 NIR 🔵 BLUE | 🛑 RED | 🔾 NIR AVAILABLE ONLY: RED/NIR AVAILABLE ONLY: BLUE/RED FACEMASK 104 PAINBUSTER 90 • Ergonomic design for application over closed eyelids. • 5"x8" neoprene flexibility for ease of application Comfortably secured by a Velcro strap. anywhere on the body. • 10.25" x 5" neoprene flexibility for ease of application · Comfortably secured by a Velcro strap. anywhere on the skin. 90 diodes of the polychromatic blend of red and infrared photon energy. • 104 diodes of the polychromatic blend of blue & red photon energy. 40 Red diodes (630 nm) in each pad. 52 Blue diodes (470 nm). • 50 NIR diodes (850 nm) in each pad. • 52 Red diodes (630 nm). 🛑 RED | 🔵 NIR 🔵 BLUE | 🛑 RED

Inlight Therapy, Inc. 2703 Broadbent Parkway NE, Suite H, Albquerque, NM 87107 • INLIGHTTHERAPYINC.COM

ILT-PRODUCTSET-521

©2021 Inlight Therapy, Inc.